

Heated Bio-Containment

Catalog # 36357

The Bio Containment module allows the housing of NHP's in a microbiological environment that can operate at a BSL-3 level independent from surrounding units or facility air handling equipment.

- All doors and openings have mechanical devices for seal compression.
- Room air is drawn into the cage through a HEPA filter on the front door.
- Removable thermometer for monitoring cage temperature.
- An Air Warming Unit can be attached to the HEPA filter to allow control of the air temperature in each cage.
- A flow indicator shows the fan is operating and the status of the filter.
- A removable fan module filters cage air before being exhausted back into the room.
- A full opening hinged cage door and sliding floor grate facilitate the entry and removal of sedated animals.

Cage Specifications:

Each Cage is 6.5 Ft.2
All Stainless Steel and Lexan Material
Suitable for 250F Autoclave

