

Transport Cart

Catalog # 35617

Constructed of Stainless Steel.

A handle is supplied to easily maneuver the cart when the bin is full.

A sturdy Rubbermaid tub holds anesthetized animals securely, and a HEPA filter is installed on the lid for airflow.

A shelf on the bottom is sloped for drainage.

The HEPA filter is locked into a hinged lid. The filter can be easily replaced and cleaned.

Cart sits on 4" Poly stem casters.

A Stainless Steel box can be made in place of the Rubbermaid tub if desired.

