

Bio-Containment Rack

Catalog # 35552

Allows the housing of NHP's in a microbiological environment that can operate at a BSL-2/3 level independent from surrounding units or facility air handling equipment.

Cage Specifications:

Each cage is 6.5 sq. ft.
Suitable for 250F autoclave

All doors and openings have mechanical devices for seal compression.

Room air is drawn into the cage through a HEPA filter on the front door.



A flow indicator shows the fan is operating and the filter's status.
A removable fan module filters cage air before being exhausted back into the room.
A full opening hinged cage door and sliding floor grate facilitate the entry and removal of sedated animals.

